

Photoluminescence mapping of $\beta\text{-Ga}_2\text{O}_3$

Cassi Remple,¹ Jesse Huso,² Ben Dutton,¹ John McCloy,¹ and Matt McCluskey^{1,2,*}

¹Washington State University, Pullman, WA 99164

²Klar Scientific, Inc., Pullman, WA 99163



WASHINGTON STATE
UNIVERSITY

Klar
Scientific

Abstract

Photoluminescence (PL) spectra of $\beta\text{-Ga}_2\text{O}_3$ include the intrinsic UV band and red emission due to Cr^{3+} impurities.¹ PL mapping with excitation wavelengths of 355 nm and 266 nm reveals the spatial distribution of these features with micron resolution. The PL maps reveal (1) striations in the red emission intensity, (2) bright UV emitters on the surface, and (3) damage induced by high-intensity laser pulses.

Method

- Klar Mini Pro PL microscope with XYZ scanning
- Ocean Optics spectrometers
- System can cover deep-UV to near-IR (266 to 2500 nm)





- T = 80 – 400 K
- This work: RT

Cr^{3+} red emission²

- (100) $\beta\text{-Ga}_2\text{O}_3\text{:Fe}$ crystal, Czochralski grown (Synoptics)
- Unintentional Cr impurities
- 355 nm excitation
- Intensity of Cr^{3+} peak (690 nm) shows variations



Striations due to variations in [Cr] or the Fermi level



(a) PL spectrum of $\beta\text{-Ga}_2\text{O}_3$ crystal. (b) PL map of the 690-nm peak intensity.

Bright UV emitters³


Bright emission at 3.27 eV localized at "pits" on the (010) surface.




PL map of the 3.27-eV peak intensity. The sample was edge-defined film-fed grown (Tamura) and hydrogenated at 900°C. As shown above, the emission occurs at surface pits.

Laser damage

- MBE-grown (010) $\beta\text{-Ga}_2\text{O}_3$ (Speck & Krishnamoorthy, UCSB), homoepitaxy
- Irradiated with 1064 nm, 3 ns pulses
- Damage threshold = 10 J/cm²
- PL excitation 266 nm



UV band emission from $\beta\text{-Ga}_2\text{O}_3$ epilayer.



PL map of the UV-band intensity for a damaged $\beta\text{-Ga}_2\text{O}_3$ epilayer. Circular "conchoidal" fracture and a cracked structure are shown.

*mattmcc@wsu.edu

Bright emitters: DOE Grant DE-FG02-07ER46386

Cr^{3+} impurities: AFOSR Grant FA9550-21-1-0078

Laser damage: LLNL Contract DE-AC52-07NA27344, LDRD Project 22-SI-003

1. M.D. McCluskey, "Point defects in Ga_2O_3 ," *J. Appl. Phys.* **127**, 101101 (2020).
2. C. Remple, J. Huso, and M.D. McCluskey, "Photoluminescence and Raman mapping of $\beta\text{-Ga}_2\text{O}_3$," *AIP Advances* **11**, 105006 (2021).
3. J. Huso, M.D. McCluskey, Y. Yu, Md. M. Islam, and F. Selim, "Localized UV emitters on the surface of Ga_2O_3 ," *Scientific Reports* **10**, 21022 (2020).